# Track® A

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PETG: copolyester

PE, insulating foil: polyethylene

Harmless to health, BfR-listed thermoplasts with tested biocompatibility.

Insoluble in water, inactive, harmless to ground water.

**Disposal/Recycling:**

Household waste (combustable, low in pollutants and high in energy)/general recycling for plastic.

**Applications:**

Invisible retainers, retention splints, active splints, occlusal splints, bruxism splints, stabilization splints, Drum splints, miniplast splints, dressing plates, moulded appliances.

**Material properties:**

Viscous hard, transparent foil. Extremely impact-resistant. Very mouldable. Bonds to acrylate.

**Finishing:**

High polish can be achieved by using plastic polishing compound. The spacer foil will be taken off only after thermoforming and finishing. Before bonding acrylic to Track A you should wet bonding area with a little bit of monomer (see brochure „Thermoforming technique“).

**Plastification:**

(for machines without heating time control) Place the foil in the devices in such a way that the spacer/insulating foil is pointing towards the model. Test softness of foil with an instrument. If permanent instrument marks show, thermoform.

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