

Aestetik-Line® Brackets

Füllstoff-verstärkte Kunststoffbrackets/Reinforced with special fillers

CE 0297

Roth-System

CE 0297

Slot .018"

MAXILLARY

	792-0503	792-0503	*792-0303	792-0203	792-0103	791-0103	791-0203	*791-0303	791-0503	791-0503
	-	-	-	-	-	-	-	-	-	-
	-	-	-	-	-	-	-	-	-	-
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angulation	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
In/Out	0,9	0,9	0,9	1,5	1,0	1,0	1,5	0,9	0,9	0,9
Rotation	2° distal	2° distal	4° mesial					4° mesial	2° distal	2° distal
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angulation	-1°	-1°	+7°	0°	0°	0°	0°	+7°	-1°	-1°
In/Out	0,7	0,7	0,9	1,5	1,5	1,5	1,5	0,9	0,7	0,7
Rotation	4° distal	4° distal	2° mesial					2° mesial	4° distal	4° distal
	792-1603	792-1503	*792-1403	791-1203	791-1203	791-1203	791-1203	*791-1403	791-1503	791-1603
	-	-	-	-	-	-	-	-	-	-

MANDIBULAR

Cases		1		5		10		20		50
Brackets		16		80		160		320		800
4-4		706-0343		706-0344		706-0345		706-0346		706-0347
Cases		1		5		10		20		50
Brackets		20		100		200		400		1000
5-5		706-0421		706-0422		706-0423		706-0424		706-0425

Slot .022"

MAXILLARY

CE 0297

	792-0501	792-0501	*792-0301	792-0201	792-0101	791-0101	791-0201	*791-0301	791-0501	791-0501
	-	-	-	-	-	-	-	-	-	-
	-	-	-	-	-	-	-	-	-	-
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angulation	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
In/Out	0,9	0,9	0,9	1,5	1,0	1,0	1,5	0,9	0,9	0,9
Rotation	2° distal	2° distal	4° mesial					4° mesial	2° distal	2° distal
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angulation	-1°	-1°	+7°	0°	0°	0°	0°	+7°	-1°	-1°
In/Out	0,7	0,7	0,9	1,5	1,5	1,5	1,5	0,9	0,7	0,7
Rotation	4° distal	4° distal	2° mesial					2° mesial	4° distal	4° distal
	792-1601	792-1501	*792-1401	791-1201	791-1201	791-1201	791-1201	*791-1401	791-1501	791-1601
	-	-	-	-	-	-	-	-	-	-

MANDIBULAR

Cases		1		5		10		20		50
Brackets		16		80		160		320		800
4-4		706-0348		706-0349		706-0350		706-0351		706-0352
Cases		1		5		10		20		50
Brackets		20		100		200		400		1000
5-5		706-0426		706-0427		706-0428		706-0429		706-0430

* Vertical-slot Ø .016" = 0,40 mm

Aestetik-Line® Brackets

Füllstoff-verstärkte Kunststoffbrackets/Reinforced with special fillers

Standard-Ricketts-System

CE 0297

Slot .018"		MAXILLARY									
		*795-0503	*795-0503	*797-0303	*797-0203	796-0103	796-0103	*796-0203	*796-0303	*795-0503	*795-0503
		-	-	-	-	-	-	-	-	-	-
		5	4	3	2	1	1	2	3	4	5
		Torque 0°	0°	+7°	+14°	+22°	+22°	+14°	+7°	0°	0°
		Angulation 0°	0°	+5°	+8°	0°	0°	+8°	+5°	0°	0°
		In/Out 1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0
		Rotation									
		Torque 0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
		Angulation 0°	0°	+5°	0°	0°	0°	0°	+5°	0°	0°
		In/Out 1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0
		Rotation									
		*795-0503	*795-0503	*796-0303	795-1203	795-1203	795-1203	795-1203	*797-0303	*795-0503	*795-0503
					-	-	-	-			
Slot .018"		MANDIBULAR									
Cases		1		5		10		20		50	
Brackets		20		100		200		400		1000	
		706-0466		706-0467		706-0468		706-0469		706-0470	

* Vertical-slot Ø .016" = 0,40 mm

Standard-Edgewise System

CE 0297

Slot .018"		MAXILLARY									
		*795-0503	*795-0503	*795-0503	795-0203	795-0103	795-0103	795-0203	*795-0503	*795-0503	*795-0503
		-	-	-	-	-	-	-	-	-	-
		5	4	3	2	1	1	2	3	4	5
		5	4	3	2	1	1	2	3	4	5
		*795-0503	*795-0503	*795-0503	795-1203	795-1203	795-1203	795-1203	*795-0503	*795-0503	*795-0503
		-	-	-	-	-	-	-	-	-	-
Slot .018"		MANDIBULAR									
Cases		1		5		10		20		50	
Brackets		16		80		160		320		800	
	4-4	706-0431		706-0432		706-0433		706-0434		706-0435	
Cases		1		5		10		20		50	
Brackets		20		100		200		400		1000	
	5-5	706-0441		706-0442		706-0443		706-0444		706-0445	