


Hybrid-Retraktor*
Hybrid retractor*

Info
N° 171




Hybrid-Retraktor: Schnelle Retraction von Eckzähnen ohne Kippung und ohne Rotation mittels superelastischer Federn aus Titanol.

SE-Federkraft: ca. 100-150 g
Federweg: 5 mm

Hybrid Retractor: Fast retraction of canines through super-elastic springs without unwanted tipping and rotation.

Retraction force: app. 100-150 g
Spring load: 5 mm

Order No.	Inhalt contents	Beschreibung description
CE 0297		
307-1020		2 Hybrid-Retraktor Hybrid retractor

Oberkiefer-/ Maxillary-kits

CE 0297

Order Nr.	Zahn/tooth	Inhalt/content
.018"	<u>3</u> <u>3</u>	2 Klebebrackets / Brackets for bonding horizontal slot .018", vertical slot .018" x .018" + 2 Hybrid retractors
.022"	<u>3</u> <u>3</u>	2 Klebebrackets / Brackets for bonding horizontal slot .022", vertical slot .018" x .018" + 2 Hybrid retractors

Unterkiefer-/ Mandibular-kits

CE 0297

Order Nr.	Zahn/tooth	Inhalt/content
.018"	<u>3</u> <u>3</u>	2 Klebebrackets / Brackets for bonding horizontal slot .018", vertical slot .018" x .018" + 2 Hybrid retractors
.022"	<u>3</u> <u>3</u>	2 Klebebrackets / Brackets for bonding horizontal slot .022", vertical slot .018" x .018" + 2 Hybrid retractors

*nach/acc. to Prof. Dr. F. G. Sander