



































Ligaturen-Elastics








Ligature elastics

Info
N° 165









Order No. universal profile CE		Farbe colour	Stück pieces	Beschreibung description
659-0008		blau/blue	50 x 20	Happy-Elastics auf Spender Happy-Elastics on key ring 
659-0009		rot/red		
659-0010		gelb/yellow		
659-0011		grün/green		
659-0012		grau/grey		
659-0013		silber/silver		
659-0014		schwarz/black		
659-0015		rosa/pink		
659-0016		marine-blau/marine-blue		
659-0017		fuchsia/pink		
659-0018		zahnfarben/tooth-colour		
659-0019		weiss/white		
659-0020		orange		
659-0021		transparent/clear		
659-0002		transparent clear		
650-0009		transparent clear	20 x 50	Elastic-Ligaturen auf Spendernadeln Elastic ligatures on cane dispensers
653-0001		transparent clear	1000	Elastic-Ligaturen Premium Qualität lose Elastic ligatures premium quality loose

Order No. universal profile CE		Farbe colour	Stück pieces	Beschreibung description
659-0100		grau/grey	10 x 100	<p>Happy-Elastics auf Pfeilspender Happy-Elastics on arrow ties</p> 
659-0101		schwarz/black		
659-0102		silber/silver		
659-0103		rot/red		
659-0104		weiss/white		
659-0105		orange		
659-0106		rosa/pink		
659-0107		fuchsia/pink		
659-0108		blau/blue		
659-0109		marine-blau/marine-blue		
659-0110		gelb/yellow		
659-0111		grün/green		
659-0112		transparent/clear		
659-0113		zahnfarben/tooth colour		
659-0116	Happy Elastics gemischt Happy Elastics mixed		10 x 100	
659-0114	Sortimentsbox (verschiede Farben) Mixed colours in Box		18 x 50	


Elastics Ketten und Zubehör
Elastics chains and accessories













Order No. universal profile		Farbe colour	Stück pieces	Beschreibung description
655-0002 small		transparent clear	1 roll 1,5 m 1 Rolle 1,5 m	 Happy-Elastics Kette auf der Rolle Happy-Elastics spool chain
655-0001 medium		transparent clear		
655-0003 wide		transparent clear		
655-0010 extra wide		transparent clear		
655-0100			10 x 22	Silicon Ligaturen Silicone ligatures
655-0015			1 Rolle / 1 Roll (1,5 m)	Lingual-Elastikkette – flacher als normale Ketten Lingual elastic chain – flatter than standard-chains

Order No.		Inhalt contents	Beschreibung description
659-0001		100	Elastic-Rotator
CE			

Order No.		Inhalt contents	Beschreibung description
657-0002 CE		1000	Separiergummiringe Ø 2,1 Separators Ø 2,1
658-0002 CE 0297		25 feet spool 7,5 m Rolle	Plasto-Flex - external Ø 0,63 mm = .025" Plasto-Flex - Außen Ø 0,63 mm = .025"
501-0869 CE		1	Ligaturen-Klemme Ligature-forceps 
501-0838 CE		1	Ligaturen-Setzer Twirl-on-Instrument 
652-0002 CE		10	Elastic-Setzer zum Setzen und Entfernen von Elastics Orthodontic wrench to attach and remove elastics 

FOR Elastics-Fruit-Line
intra-oral



LATEX Elastics	NON-LATEX Elastics		Inhalt contents	Ring inner dia. mm	inch	Beschreibung description
CE	CE					
650-1006	650-4006		10 x 100	3,2	1/8	schwach light 60 g = 2 OZ
650-1001	650-4001		10 x 100	4,8	3/16	
650-1002	650-4002		10 x 100	6,3	1/4	
650-1003	650-4003		10 x 100	7,9	5/16	
650-1016	650-4016		10 x 100	3,2	1/8	mittel medium 115 g = 4 OZ
650-1011	650-4011		10 x 100	4,8	3/16	
650-1012	650-4012		10 x 100	6,3	1/4	
650-1013	650-4013		10 x 100	7,9	5/16	
650-1026	650-4026		10 x 100	3,2	1/8	heavy stark 170 g = 6 OZ
650-1021	650-4021		10 x 100	4,8	3/16	
650-1022	650-4022		10 x 100	6,3	1/4	
650-1023	650-4023		10 x 100	7,9	5/16	

For more details see www.forestadent.com