

# Aestetik-Line® Brackets

Füllstoffverstärkte Kunststoffbrackets/Reinforced with special fillers

## Slot .018"

### Roth-System



		MAXILLARY										CE 0297
		792-0503	792-0503	*792-0303	792-0203	792-0103	791-0103	791-0203	*791-0303	791-0503	791-0503	
		-	-	-	-	-	-	-	-	-	-	
		Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
		Angulation	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
		In/Out	0,9	0,9	0,9	1,5	1,0	1,0	1,5	0,9	0,9	0,9
		Rotation	2° distal	2° distal	4° mesial					4° mesial	2° distal	2° distal
		Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
		Angulation	-1°	-1°	+7°	0°	0°	0°	0°	+7°	-1°	-1°
		In/Out	0,7	0,7	0,9	1,5	1,5	1,5	1,5	0,9	0,7	0,7
		Rotation	4° distal	4° distal	2° mesial					2° mesial	4° distal	4° distal
		792-1603	792-1503	*792-1403	791-1203	791-1203	791-1203	791-1203	*791-1403	791-1503	791-1603	
		-	-	-	-	-	-	-	-	-	-	

### MANDIBULAR

<b>Fälle / Cases</b>	1	5	10
Brackets	16	80	160
4-4	706-0343	706-0344	706-0345
<b>Fälle / Cases</b>	1	5	10
Brackets	20	100	200
5-5	706-0421	706-0422	706-0423

## Slot .022"



		MAXILLARY										CE 0297
		792-0501	792-0501	*792-0301	792-0201	792-0101	791-0101	791-0201	*791-0301	791-0501	791-0501	
		-	-	-	-	-	-	-	-	-	-	
		Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
		Angulation	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
		In/Out	0,9	0,9	0,9	1,5	1,0	1,0	1,5	0,9	0,9	0,9
		Rotation	2° distal	2° distal	4° mesial					4° mesial	2° distal	2° distal
		Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
		Angulation	-1°	-1°	+7°	0°	0°	0°	0°	+7°	-1°	-1°
		In/Out	0,7	0,7	0,9	1,5	1,5	1,5	1,5	0,9	0,7	0,7
		Rotation	4° distal	4° distal	2° mesial					2° mesial	4° distal	4° distal
		792-1601	792-1501	*792-1401	791-1201	791-1201	791-1201	791-1201	*791-1401	791-1501	791-1601	
		-	-	-	-	-	-	-	-	-	-	

### MANDIBULAR

<b>Fälle / Cases</b>	1	5	10
Brackets	16	80	160
4-4	706-0348	706-0349	706-0350
<b>Fälle / Cases</b>	1	5	10
Brackets	20	100	200
5-5	706-0426	706-0427	706-0428

\* Vertical-slot Ø .016" = 0,40 mm

# Aestetik-Line® Brackets

Füllstoffverstärkte Kunststoffbrackets/Reinforced with special fillers

**Slot .018"**

## Standard-Ricketts-System

CE 0297



		MAXILLARY									
		*795-0503	*795-0503	*797-0303	*797-0203	796-0103	796-0103	*796-0203	*796-0303	*795-0503	*795-0503
		-	-	-	-	-	-	-	-	-	-
Torque		0°	0°	+7°	+14°	+22°	+22°	+14°	+7°	0°	0°
Angulation		0°	0°	+5°	+8°	0°	0°	+8°	+5°	0°	0°
In/Out		1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0
Rotation											
Torque		0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Angulation		0°	0°	+5°	0°	0°	0°	0°	+5°	0°	0°
In/Out		1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0	1,0
Rotation											
		*795-0503	*795-0503	*796-0303	795-1203	795-1203	795-1203	795-1203	*797-0303	*795-0503	*795-0503
		-	-	-	-	-	-	-	-	-	-

MANDIBULAR

Fälle / Cases	1	5	10
Brackets	20	100	200
	706-0466	706-0467	706-0468

\* Vertical-slot Ø .016" = 0,40 mm

**Slot .018"**

## Standard-Edgewise System

CE 0297



		MAXILLARY									
		*795-0503	*795-0503	*795-0503	795-0203	795-0103	795-0103	795-0203	*795-0503	*795-0503	*795-0503
		-	-	-	-	-	-	-	-	-	-
		*795-0503	*795-0503	*795-0503	795-1203	795-1203	795-1203	795-1203	*795-0503	*795-0503	*795-0503
		-	-	-	-	-	-	-	-	-	-

MANDIBULAR

Fälle / Cases	1	5	10	
Brackets	16	80	160	
	4-4	706-0431	706-0432	706-0433
Fälle / Cases	1	5	10	
Brackets	20	100	200	
	5-5	706-0441	706-0442	706-0443