






# Ni-Free Products





Bukkalröhrchen zum Kleben – Ni-Free

Buccal tubes for direct bonding – Ni-Free

## OK Maxillary

	Zahn tooth	Torque	Angulation	distaler Offset distal offset		 Pad sandgestraht Sandblasted pad auf Pad (10 Stück/10 pcs.) rechts/right   links/left
 convertible		-14°	0°	14°	.018" x .025"	741M0704   740M0704
					.022" x .028"	741M0702   740M0702
		-14°	0°	14°	.018" x .025"	741M0804   740M0804
					.022" x .028"	741M0802   740M0802

## UK Mandibular

 convertible		-30°	0°	4°	.018" x .025"	741M1704   740M1704
					.022" x .028"	741M1702   740M1702
		-30°	0°	4°	.018" x .025"	741M1804   740M1804
					.022" x .028"	741M1802   740M1802

\*Nickel free