


Mechanical retention is better!



Aspire Gold Ceramic Brackets


NEW!
Now also available in McLaughlin-
Bennett-Trevisi System*




- **Handles like a metal bracket...**
The ASPIRE GOLD ceramic bracket convinces through its natural look and low profile. Undercuts in the bracket base give excellent retention, similar to metal brackets with a mesh base. There is no need for the use of special bonding enhancements. The bracket can be bonded with standard light or chemical cure bonding agents.
- **And debonds like a metal bracket...**
no need for special debonding instruments or a special technique.




Aspire Gold Ceramic Brackets



Compared to an ordinary ceramic-bracket the gold slot of the Aspire Gold bracket helps to reduce the friction to a minimum.



Surface of a BioStarter® NiTi-archwire in River Finish Gold



Surface of a Standard NiTi-archwire

To reduce the friction even further you can use the river finish gold range of NiTi-wires. The surface of the wires is specially treated with a plasma coating.



The river finish gold coating is standard in the Forestadent range of BioStarter, Low Force and Triple Force archwires and helps to enhance the aesthetic appearance of the appliance.

Slot .018"

Roth-System

CE 0297

Cases	5
Brackets 3-3	60
706-0731	
3 706-0732	
-	-

Cases	1	5	10	20	50
Brackets 5-5	20	100	200	400	1000
706-0700	706-0701	706-0702	706-0703	706-0704	
3 706-0710	706-0711	706-0712	706-0713	706-0714	
3-5 706-0720	706-0721	706-0722	706-0723	706-0724	

Aspire Gold Brackets are available per piece

Slot .022"

Roth-System

CE 0297

Cases	5
Brackets 3-3	60
706-0733	
3 706-0734	
-	-

Cases	1	5	10	20	50
Brackets 5-5	20	100	200	400	1000
706-0705	706-0706	706-0707	706-0708	706-0709	
3 706-0715	706-0716	706-0717	706-0718	706-0719	
3-5 706-0725	706-0726	706-0727	706-0728	706-0729	

Aspire Gold Brackets are available per piece

Slot .022"

McLaughlin-Bennett-Trevisi-System*

NEW!

CE 0297

Cases	5
Brackets 3-3	60
706-0750	
3 706-0751	
-	-

Cases	1	5	10	20	50
Brackets 5-5	20	100	200	400	1000
706-0735	706-0736	706-0737	706-0738	706-0739	
3 706-0740	706-0741	706-0742	706-0743	706-0744	
3-5 706-0745	706-0746	706-0747	706-0748	706-0749	

*The FORESTADENT version of these prescriptions are not claimed to be a duplication of any other, nor does Forestadent imply that it is endorsed in any way by the doctor.

Aspire Gold Brackets are available per piece