

Slim Line Expanders

Compact & Multi-Functional

Forestadent has developed a versatile and compact expansion screw with 2 arms. This expansion screw offers an outstanding solution for both upper and lower appliances. The compact design facilitates while allowing earlier expansion of the maxillary arch.



Quality System



Certified by DOS acc.to
DIN EN ISO 9001 : 2000
DIN EN ISO 13485 : 2003
Annex II of Directive
93/42/EEC
Reg. No. 55387



Part Numbers:

8mm #198-1223

12mm #198-1623

Blue Expansion Keys #100-0009

Lighted Expansion Keys #100-10010

- Both expanders provide uprighting and torquing of molars, expansion of the maxillary arch in the posterior segments, distalization and the possibility of mesializing of the posterior region.
- Extremely strong stainless steel material provides guaranteed results. The spindles produce less friction during treatment.
- The body of the screw is extremely small and is anatomically designed to fit extremely small palates.
- The screws are available in two unique designs, either 8mm or 12mm total expansion. The body of the screws are permanently identified by a laser marked arrow.
- Besides providing bilateral expansion, the new Slim Line is easy to use for both patient and clinician.
- One complete turn or 360° equals 0.8mm of expansion.

Contact Your Local Representative or
Call Toll Free at **1-800-721-4940**
The Next Level in Orthodontics

