



Strength, Durability, Reliability.

Just a few important features of our newly designed expansion screws



Palatal Split Screw

167-1539
Pre-bent arms
12 mm
15 mm wide



Palatal Split Screw

167-1326
Pre-bent arms
10 mm
13 mm wide



Palatal Split Screw

167-1323
Straight arms
10 mm
13 mm wide



Palatal Split Screw

169-1323
For acrylic use
10 mm
13 mm wide

Lower Lingual Expander



169-1323S
10 mm
13 mm wide



Lower Lingual Expander

168-1223
7 mm
11 mm wide

We have taken rapid palatal expansion to the next level with these two totally new designs. By using the highest quality non-corrosion stainless steel available and pre-bending the extension arms you will have the most durable and reliable screw on the market. We also have inserted the extension arms inside the body of the screw for strength during expansion.

Forestadent's long history of high-quality expansion screws guarantees consistency in your practice.



Contact Your Local Representative or Call Toll Free at **1-800-721-4940** • The Next Level In Orthodontics